HANDHABUNGSHINWEIS
INSTRUCTIONS FOR USE
ISTRUZIONI PER L'USO
ISTRUCCIONES DE USO

CLASSIC

Bipol. coagulation forceps
Pincés à coag. bipolaire
Pinza per coag. bipolare
Pinza de coag. bipolar

85500000 - 85588013

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English Bipolar coagulation forceps

Attention
Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Intended Use
The CLASSIC bipolar coagulation forceps is intended for use in minimally invasive and especially laparoscopic surgical procedures. The device is to be passed through a 5.5 mm laparoscopic cannula. The CLASSIC bipolar coagulation forceps is designed to manipulate, grasp or cut selected tissue. It must be connected to the bipolar output of an electrosurgical generator. As indicated bipolar coagulation current may be selectively applied to the tissue.

Maximum output voltage of the generator. U_max: 500 V_p

Appropriate connecting cables/adapters:
Bissinger bipolar cable REF 801 0010x.
Bissinger adapter REF 85500400, 85500401

Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.
Caution: In particular when using the scissors, parenchymal tissue may deflagrate.

Contraindications
Incidents which have been reported in connection with the use of bipolar systems:
Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.
Fire in connection with surgical drapes and other inflammable materials.
Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
Explosions caused by sparks in the proximity of inflammable gases.
Perforation of organs. Sudden severe bleedings.

Use and safety instructions
The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- The instrument may not be laid down on the patient.
- Coagulation should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.
Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

Assembly and Operation
Follow the pictogram PIC-855 for assembly and disassembly of the instrument.
Once correctly assembled and connected to power, the device may be used in either the right or the left hand.
To close jaws: compress (grip) handle.
To open jaws: release (grip) handle.
Coagulation current is activated by a footswitch that is part of your electrosurgical unit.

Reprocessing
Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.
Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport
Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

- Machine reprocessing
Cleaning
Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.
1. Pre-rinse for 1 min. with cold water
2. Discharging
3. Pre-rinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).
10. Discharging

Disinfection
Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying
Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.
If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing
Ultrasound pre-cleaning
1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasound must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

Cleaning
Prepare a cleaning bath according to the manufacturer's instructions.

- 1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
4. The following step only applies to channels and the insides of tubes: Push the brush into and out of the tubes at least six times. Rinse the tubes with DI water. Repeat the procedure.
5. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection
Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying
Manual drying is carried out using a lintfree cloth and, in particular, for drying cavities and channels, sterile compressed air.

Functional test and packaging
Perform visual inspection for cleanliness; if required, perform an assembly and functional test according to the operating instructions.
If necessary, repeat the reprocessing process until the instrument is optically clean.
Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation
Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

Storage
Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

Attention: Defect products must pass the complete reprocessing process before being returned for repair.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):

- Neodisher FA by Dr. Weigert (alkaline)
- Endozime by Ruhof (enzymatic)

Cleaning agents (manual cleaning):

- Enzol Enzym, detergent by Johnson&Johnson
- Disinfectants (manual disinfection):

- Cidex OPA, Johnson&Johnson
- Neutrising agent:

- Neodisher Z by Dr. Weigert
- Cleaning and disinfection device:

- Miele G 7736 CD
- Miele insert module E 327-06
- Miele MIS module E 450

For details, see report.

- SMP GmbH # 01707011901 (machine cleaning)
- MDS GmbH # 135196-10 (man. cleaning/disinfection)
- Nelson Labs # 200432706-02 (sterilisation)
- MDS GmbH Testbericht 084183-10 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexpert way.

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)
DQS Medizinprodukte GmbH