# bissinger

HANDHABUNGSHINWEIS  
INSTRUCTIONS FOR USE  
INSTRUCTIONS DE SERVICE  
ISTRUZIONI PER L'USO  
INSTRUCCIONES DE USO

BiTech™

Bipolare Scheren  
Bipolar Scissors  
Ciseaux bipolaires  
Forbici bipolari  
Tijeras bipolares

REF

85713150 – 85828120

CE 0297

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## ENGLISH BiTech™ Bipolar Scissors

### Caution

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments or may pose a hazard to patients and users.

### Intended Use

The Bissinger bipolar scissors are designed to cut, dissect and coagulate selected tissue.

They must be connected to the bipolar output of an electro-surgical generator using a suitable bipolar cable and must only be used with parameters for bipolar coagulation.  
Do not exceed a maximum output of 300 V<sub>p</sub> of your generator.

### Suitable connecting cables:

Bissinger Bipolar Cables REF 80100081 - 80100088

**Attention:** Instruments for electro-surgery should be used only by persons who have been specially trained in the use of such instruments.

### Contraindications

**Adverse events reported while using bipolar electro-surgical devices include:**

- Inadvertent activation with resultant tissue damage at the wrong site and/or equipment damage.
- Fires involving surgical drapes and other combustible materials have been reported.
- Alternate current pathways resulting in burns where the patient or physician or assistant is in contact with exposed metal.
- Explosions caused by electro-surgical sparking in a flammable gas environment (i.e. explosive anesthetic gases)
- Organ perforation. Sudden massive hemorrhage.

### Use and safety instructions

- All instruments have to be completely cleaned, disinfected and sterilised before initial use and any other use.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular, areas such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.
- **Bipolar scissors include high-quality ceramic parts which have to be treated with particular care and protected against breakage.**
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- The instrument may not be laid down on the patient.
- Coagulation should only be performed if the contact surfaces are visible. Do not touch any other metallic instruments during coagulation.

### Cautions & Warnings

**US Federal Law restricts this device to sale by or on the order of a physician only.**

Special safety precautions should be observed when using electro-surgical instruments.  
Electro-surgical instruments can pose a significant shock, burn or explosion hazard if used improperly, incorrectly or carelessly.

Avoid touching or grounding electro-surgical instruments to non-insulated instruments, scopes, trocar sleeves, etc. All persons using such devices should be knowledgeable in the use and handling of laparoscopic instruments, coagulation equipment, their accessories and other related equipment.  
Test all instruments, accessories and equipment prior to each use. Do not use in presence of flammable liquids or anaesthetics.

Electro-surgical generators used with these devices are designed to cause destruction of tissue and are inherently dangerous if operated improperly. Follow all safety precautions and instructions supplied by the manufacturer of the electro-surgical generator.

The electrode tip must always be in full view before activating power. Apply power only when electrode tip is in full contact with the tissue selected for coagulation. Electrode tip must not come in contact with other metal instruments during use.  
Failure to observe these cautions and contraindications may result in injury, malfunction or other unanticipated occurrences or events for the operator, staff and/or the patient.

### Reprocessing



Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.

Instruments for electro-surgery are by nature subject to increased wear depending on the type and time of use.

### Preparation and transport

Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C) as this may result in the fixation of residues and could reduce the cleaning success.

Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container to avoid any damage to the instruments and any contamination of the environment.  
Handle with extreme care, do not throw or drop !

### Machine reprocessing

#### Manual reprocessing

1. Insert the Instrument for 5 minutes in cold tap water.
2. Clean the Instruments under cold tap water with a soft brush until all visible dirt is removed.

#### Cleaning

Place the instrument with the scissor blades open in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Pre-rinse for 1 min. with cold water
2. Discharging
3. Pre-rinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).
10. Discharging

#### Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

#### Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.  
If necessary, manual drying may additionally be carried out using a lint free cloth. Dry cavities of the instruments by blowing with sterile compressed air.

### Manual reprocessing

#### Ultrasonic Pre-Cleaning

1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasonic must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

#### Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.

3. Clean the instrument in the bath manually using a soft brush. All surfaces have to be brushed several times.
4. Rinse the products thoroughly with running tap water to remove the cleaning agents without residue.

#### Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

#### Drying

Manual drying is carried out using a lintfree cloth and, in particular, for drying cavities and channels, sterile compressed air.

#### Functional test and packaging

Perform visual inspection for cleanliness; if required, perform an assembly and functional test according to the operating instructions. If necessary, repeat the reprocessing process until the instrument is optically clean. Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

#### Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.

- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

#### Storage

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

#### Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

**Attention:** Defect products must pass the complete reprocessing process before being returned for repair.

#### Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

#### Cleaning agents (for machine use):

Neodisher FA by Dr. Weigert (alkaline)  
Endozime by Ruhof (enzymatic)

#### Cleaning agents (manual cleaning):

Enzol Enzym, detergent by Johnson&Johnson

#### Disinfectants (manual disinfection):

Cidex OPA, Johnson&Johnson

#### Neutralising agent:

Neodisher Z by Dr. Weigert

#### Cleaning and disinfection device:

Miele G 7736 CD  
Miele insert module E 327-06  
Miele MIS module E 450

For details, see report.

SMP GmbH # 01707011901 (machine cleaning)  
MDS GmbH # 135196-10 (man. cleaning/disinfection)  
Nelson Labs # 200432706-02 (sterilisation)  
MDS GmbH Testbericht 084183-10 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

#### Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

#### Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexperienced way.

#### Explanation of symbols

	LOT-Number
	Caution: Non-sterile product
	Reference number
	Attention !
	Refer to instructions for use
	CE-Mark and registration number of the Notified Body DQS Medizinprodukte GmbH August-Schanz-Straße 21 60433 Frankfurt, Germany
	Manufacturer Production date

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)