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English Bipolare Resectoscope

Attention

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Intended Use

The PLASMALOOP was developed in particular for the transurethral resection (TUR). The instrument used to ablate in layers and coagulate from tissue of the bladder or prostate gland. Using of suitable electrode (REF 867000xx) tissue can be cut or coagulate.
Using of a Resectoscope sheath, continuously rinsing sheath (CH24/26; CH24/27) or single rinsing sheath (CH24), the PLASMALOOP is introduced. The cutting or coagulating result from electric energy, generated by HF generators for the electric surgery. Using of a suitable bipolar cable it is connected to the bipolar exit of a HF generator. If indicates, bipolar coagulations or cut current can be used straight.

Maximum output voltage of the generator. U_{max}:

800 V_p

Suitable connector cable:

Bissinger Bipolar cable REF 801 00xxx.

Suitable accessories:

Resectoscope shaft REF 867100xx
Cystoscope REF 30-0xxx-xx

Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.

Contraindications

Incidents which have been reported in connection with the use of bipolar systems:

Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.
Fire in connection with surgical drapes and other inflammable materials.

Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.

Explosions caused by sparks in the proximity of inflammable gases.
Perforation of organs. Sudden severe bleedings.

Use and safety instructions

The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.

- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts and insulations have to be checked carefully.
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- The instrument may not be laid down on the patient.
- Resection may only be used in an electrically leading-capable physiological salt solution.
- Coagulation and cut should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics, shafts or things like that during use.

Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

Assembly and Operation

Follow the pictogram PIC-867 for assembly and disassembly of the instrument.
Once correctly assembled and connected to power, the device may be used in either the right or the left hand. With operate of the handle the electrode is moved.
Coagulation current is activated by a footswitch that is part of your electrosurgical unit.

Reprocessing

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.

Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport

Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

Complex instruments must be taken apart for cleaning and disinfection in accordance with pictogram.

Manual Pre-Cleaning

1. Immerse the instrument in cold water for 5 minutes. If possible disassemble the instruments and clean with a soft brush under cold water until all visible impurities are removed. In cavities, holes and threads flush with a water jet pistol at least for 10 sec. (pulsed process).
2. Place instrument in an ultrasonic bath with 0.5 % alkaline (enzymatic) cleaning detergent. Ultrasound must be applied for 15 minutes at 40°/104°F
3. Remove the instrument and rinse completely with cold tap water to remove the cleaning detergent.

Machine reprocessing

Cleaning

Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Prerinse for 1 min. with cold water
2. Discharging
3. Prerinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).
10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.
If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Cannot be applied for this instrument.

Functional test and packaging

Perform visual inspection for cleanliness good condition of the materials. If necessary, repeat the reprocessing process until the instrument is optically clean.
Perform an assembly and functional test. Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

Storage

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

Attention: Defect products must pass the complete reprocessing process before being returned for repair.

Information on the validation of the reconditioning information on the validation of the reconditioning
The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):

Neodisher FA by Dr. Weigert (alkaline)
Endozime by Ruhof (enzymatic)

Cleaning agents (manual cleaning):

Enzol Enzym, detergent by Johnson&Johnson

Disinfectants (manual disinfection):

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If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexperienced way.

Explanation of symbols

LOT-Number

Caution: Non-sterile product

Reference number

Attention !

Refer to instructions for use

CE-Mark and registration number of the Notified Body
DQS Medizinprodukte GmbH
August-Schanz-Straße 21
60433 Frankfurt, Germany

Manufacturer
Production date

① Assembly
② Disassembly

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SMP GmbH # 01707011901 (Masch. Reinigung)
MDS GmbH Testbericht 084183-10 (Sterilisation)