bissinger

HANDHABUNGSHINWEIS
INSTRUCTIONS FOR USE
INSTRUCTIONS DE SERVICE
ISTRUZIONI PER L'USO
ISTRUCCIONES DE USO

Elektrodenhandgriffe
Handles for electrodes
Poignees pour electrodes
Impugnatura per elettrodi
Mangos porta electrodos

REF

89213022, 89213040,
89214040, 89215040,
89216040, 89217040,
89218040, 89214024,
89215024, 89216024,
89217024, 89218024,
89219024, 89216025

CE 0297

Günter Bissinger Medizintechnik GmbH
Hans-Theisen-Str.1
79331 Teisingen
Germany

Tel.: ++49 7641 9 14 33 0
Fax: ++49 7641 9 14 33 33
Email: info@bissinger.com
www.bissinger.com

HH-892_Elektrodenhandgriffe-2015multilingual_vC.docx
Revision C 04.05.2015/MB

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Achtung

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English Handles for electrodes

Attention

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Intended Use

The HF electrode handles are made for cutting and coagulating biological tissue. They are used to hold shaft electrodes of any known electrode shape with a shaft Ø 4 mm (± 0.1 mm) or Ø 2.4 mm (± 0.1 mm) (depending on design). The electrode can be inserted twist-proof in different positions.
The handle is connected to the monopolar output of an HF generator with a corresponding cable and must only be used with the intended parameters.
Maximum power output allowed Umax: 5 kVp

The cutting or coagulation current is set at the HF generator. It is triggered via a button at the handle (yellow button = activation of the cutting current, blue button = activation of the coagulation current or grey button = activation of the mode set at the device). Triggering is also possible via the foot switch connected to the HF generator.

Suitable monopolar electrodes include:

- for shaft Ø4 mm: Bissinger REF 89500xxx
- for shaft Ø 2.4 mm: Bissinger REF 89503xxx, 89505xxx

Suitable connecting cables:

Bissinger REF 89101071, 89101072.

Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.

Contraindications

Incidents which have been reported in connection with the use of electrosurgical systems:
Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.
Fire in connection with surgical drapes and other inflammable materials.
Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
Explosions caused by sparks in the proximity of inflammable gases.
Perforation of organs. Sudden severe bleedings.
Bipolar forceps have proved inefficient for tubular sterilisation or coagulation in the context of sterilisation and should therefore not be used for this purpose.

Use and safety instructions

The non-obsolescence of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- The instrument may not be laid down on the patient.
- Coagulation should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.
- Do not use in patients with pacers or other active implants without seeking the corresponding technical advice first.
- Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.
Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

Reprocessing

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.
Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport

Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

Machine reprocessing

Cleaning
Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

- 1. Pre-rinse for 1 min. with cold water
2. Discharging
3. Pre-rinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).
10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine. If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Ultrasonic Pre-Cleaning

- 1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasound must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

- 1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
4. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Manual drying is carried out using a lintfree cloth and, in particular, for drying cavities and channels, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness; if required, perform an assembly and functional test according to the operating instructions.
If necessary, repeat the reprocessing process until the instrument is optically clean.
Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C.
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

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Neodisher FA; Dr. Weigert (alcalino)
Endozime; Ruhof (enzimatico)

Detergente (manuale):

Enzol Enzym Detergent, Johnson&Johnson

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)
DQS Medizinprodukte GmbH